

Form B03

Scheduled Maintenance Work Order



Quantum Medical Solutions sdn bhd
transcending boundaries, transforming life

Format Ref: - QMS/TSD-022 Rev.01

Work Order No.	PWO346736	Schedule Month	November 2018
Work Order Date	01/11/2018	Completed Date	27/11/18
Clinic Name	Klinik Pergigian Lundu, Hospital Lundu	Clinic Code	SWK310
BE No.	SWK004932	Distict	KUCHING
BE Category	Micromotor	WO Assigned to	Foster Anak Engkasan
Ownership	<input checked="" type="checkbox"/> Existing Equipment	<input type="checkbox"/> Purchase	<input type="checkbox"/> New
BE Condition	<input checked="" type="checkbox"/> Active	<input type="checkbox"/> BER Proposed	
Work Order Type	<input checked="" type="checkbox"/> Preventive Maintenance (PM)	<input type="checkbox"/> Third Party Calibration (TPC)	
	<input type="checkbox"/> Routine Inspection (RI)	<input type="checkbox"/> Statutory Certification (SC)	

Reschedule Date	
-----------------	--

BE Third Party Calibration / Statutory Certification Details

Company Name		Cal / Cert Date	
Contact Number		Cal / Cert Expiry Date	

Action Taken

- to do ppm perbeck's attend.
- ppm done.

Schedule Maintenance Execution Details

SI No	QMS Engineer / Technician Name	Date	Start Time	End Time
01	FOSTER ANAK ENKASAN QUANTUM BIOMEDICAL ENGINEERING KUCHING SARAWAK.	27/11/18	1000	1020

Customer Remarks

Engineer / Technician Signature Name FOSTER ANAK ENKASAN Date QUANTUM BIOMEDICAL ENGINEERING KUCHING SARAWAK. 27/11/18	Customer Signature Name DOREENA SALIPUI Designation Juruterapi Pergigian U20 Date 27/11/18 Seal
--	--

For Internal Use

First Verification
QMS Circle Incharge

Final Verification
QMS State Incharge



KEMENTERIAN KESIHATAN MALAYSIA

MEET Planned Preventive Maintenance Checklist

Micromotor

BE CODE : DE-019

CHECKLIST NO: GL-096-000
REV.000

PART 1 ASSET DETAILS

WORK ORDER NO	▶ P06346736	BE NO	▶ SWK004932
MANUFACTURER	▶ NSK, Japan	MODEL	▶ M10
FREQUENCY	▶ 6 MONTHLY () 12 MONTHLY (✓)	PPM HOURS	▶ 1.00

PART 2 SPECIAL PRECAUTION

If there is evidence of body fluid contamination, submit the device for cleaning and decontamination before inspecting it.

Wear appropriate Personnel Protection Equipment (PPE) during work.

Wear grounded electrostatic wristband when handling PCB or electronic components.

Refer to the safety procedure for additional precautions and guidance as per manufacturer guidelines.

Make sure the test equipment used are duly calibrated.

PART 3 TEST APPARATUS

Tick (✓) where appropriate

NO	ASSET NO	DESCRIPTION	SERIAL NO	CALIBRATION DUE ON
		ELECTRICAL SAFETY ANALYZER		
		TACHOMETER		

PART 4 QUALITATIVE TASKS

Tick (✓) where appropriate

	PASS	FAIL	NA		PASS	FAIL	NA
1 Chassis - verify physical integrity, cleanliness and condition	(✓)	()	()	8 Label - verify physical integrity	(✓)	()	()
2 Mount/Fasteners - verify physical integrity	(✓)	()	()	9 Indicators/ Displays - verify proper illumination and operation	(✓)	()	()
3 AC Plug / Power Cord - verify physical integrity and proper insulation	(✓)	()	()	10 Motor - verify proper operation physical integrity	(✓)	()	()
4 Strain Relief - verify physical integrity at both ends of line cord	(✓)	()	()	11 Accessories-verify physical integrity and operations.	(✓)	()	()
5 Circuit Breaker/ Fuse - verify integrity of external circuit breaker and/or rating of external fuse	(✓)	()	()	12 Foot/Knee control-verify proper operation	(✓)	()	()
6 Fittings/ Connectors - check all fittings/connectors	(✓)	()	()				
7 Controls/Switches - verify proper operation of controls	(✓)	()	()				

PART 5 PREVENTIVE MAINTENANCE TASKS

Tick (✓) where appropriate

	DONE	NOT DONE **	NA		DONE	NOT DONE **	NA
1 Cleanliness - clean interior and exterior of the equipment	(✓)	()	()	4 Check micromotor speed	(✓)	()	()
2 Service micromotor head/chuck	(✓)	()	()	Notes: * For all Parts, NA is defined as NOT APPLICABLE			
3 Align / adjust mechanical components	(✓)	()	()	*** If you have ticked 'NOT DONE', then justify in Part 8			
				*** Choose Whichever Applicable			