

Scheduled Maintenance Work Order

CMP

Quantum Medical Solutions sdn bhd
transcending boundaries. transforming life.

Format Ref:- QMS/TSD-022 Rev.01

Work Order No	PW0410387	Scheduled Month	January 2019
Work Order Date	01/01/19	Completed Date	9/1/19
Clinic Name	KL Alor Pongsu	Clinic Code	PRK018
BE No	PRK030704	District	Kerian
BE Category	Sphygmomanometer, Aneroid	WO Assigned to	Razin
Ownership	<input checked="" type="checkbox"/> Existing Equipment	<input type="checkbox"/> MEET Equipment	
BE Condition	<input type="checkbox"/> Under Warranty	<input type="checkbox"/> BER Proposed	
Work Order Type	<input type="checkbox"/> Preventive Maintenance (PM) <input checked="" type="checkbox"/> Routine Inspection (RI)	<input type="checkbox"/> Third Party Calibration (TPC) <input type="checkbox"/> Statutory Certification (SC)	

Reschedule Date	NA
-----------------	----

BE Third Party Calibration / Statutory Certification Details

Company Name	NA	Cal / Cert Date	NA
Contact Number	NA	Cal / Cert Expiry Date	NA

Action Taken

PPM done as per checklist

Manufacturer : Welch Allyn
 Model : Shock Resistant
 S/N : 14P11183207
 Purchase Date : NA
 Kewpa No : NA

Schedule Maintenance Execution Details

SI No	QMS Engineer / Technician Name	Date	Start Time	End Time
1	Razin	9/1/19	1430	1500

Customer Remarks

Engineer / Technician Signature Name Date MUHAMMAD RAZIN BIN MOHAMAD RAZIN JUNIOR BIOMEDICAL TECHNICIAN QUANTUM MEDICAL SOLUTIONS SDN BHD 9/1/19	Customer Signature Name Designation Date Seal NURUL HASANAH BT HUSAIN JURURAH 9/1/19
--	---

For Internal Use only

First Verification
QMS Circle InchargeFinal Verification
QMS State Incharge



KEMENTERIAN KESIHATAN MALAYSIA

MEET Planned Preventive Maintenance Checklist
Sphygmomanometers, Aneroid
BE CODE : 16-156

CHECKLIST NO: CL-140-000
REV.000

PART 1 ASSET DETAILS

WORK ORDER NO ▶ PWO 410387 BE NO ▶ PRK030704
MANUFACTURER ▶ Welch Allyn MODEL ▶ Shock Resistant
FREQUENCY ▶ 12 MONTHLY (✓) PPM HOURS ▶ 0.50

PART 2 SPECIAL PRECAUTION

If there is evidence of body fluid contamination, submit the device for cleaning and decontamination before inspecting it.

Wear appropriate Personnel Protection Equipment (PPE) during work.

Wear grounded electrostatic wristband when handling PCB or electronic components.

Refer to the safety procedure for additional precautions and guidance as per manufacturer guidelines.

Make sure the test equipment used are duly calibrated.

PART 3 TEST APPARATUS

Tick (✓) where appropriate

NO	ASSET NO	DESCRIPTION	SERIAL NO	CALIBRATION DUE ON
1	NA	PRESSURE METER	NA	NA
2	TEESA0262	Pressure Gauge	002040	2/12/19

PART 4 QUALITATIVE TASKS

Tick (✓) where appropriate

	PASS	FAIL	NA
1 Chassis - Verify physical integrity, cleanliness and condition.	(/)	()	()
2 Mount/Fasteners - Verify physical integrity.	(/)	()	()
3 Check cuff & hose - Verify physical integrity and cleanliness.	(/)	()	()
4 Check Bulb - Verify physical integrity and cleanliness.	(/)	()	()
5 Check air release valve.	(/)	()	()
6 Check dial.	(/)	()	()

PART 5 PREVENTIVE MAINTENANCE TASKS

Tick (✓) where appropriate

	DONE	NOT DONE **	NA
1 Clean the Exterior/Interior.	(/)	()	()

Notes:

- * For all Parts, NA is defined as NOT APPLICABLE
- ** If you have ticked 'NOT DONE', then justify in Part 8
- *** Choose Whichever Applicable



KEMENTERIAN KESIHATAN MALAYSIA

MEET Planned Preventive Maintenance Checklist
Sphygmomanometers, Aneroid
BE CODE 16-156

CHECKLIST NO: CL-140-000
REV.000

WORK ORDER NO ▶ PWO 10387

PART 6 QUANTITATIVE TASKS

Tick (✓) where appropriate

No	Description	Units / UOM	Set Values	Measured Values	Limit/Tolerance	PASS	FAIL	NA
1	Blood pressure accuracy	mmHg	40	110	38-42	(/)	()	()
		mmHg	70	70	68-72	(/)	()	()
		mmHg	100	100	98-102	(/)	()	()
		mmHg	130	130	128-132	(/)	()	()
		mmHg	160	160	156-162	(/)	()	()
		mmHg	190	190	188-192	(/)	()	()

PART 7 ELECTRICAL SAFETY TEST

ELECTRICAL SAFETY TEST, (attach report)

Tick (✓) where appropriate

Standard use:

☐ IEC 60601

☐ IEC 61010

☐ IEC 62353

Result:

☐ PASS

☐ FAIL

☒ NA

PART 8 NOTES

* Est not required due to not electrical usage.



CORRECTIVE MAINTENANCE REQUIRED



FUNCTIONING



NOT FUNCTIONING

WORK ORDER NO ▶

NA

NEXT PPM DATE ▶

Jan 20

PPM has been performed in accordance to the checklist and the equipment is functioning to the intended purpose.

COMPLETED BY:

MUHAMMAD NAZIB BIN MOHAMMAD
JUNIOR ENGINEER
QUANTUM MEDICAL SOLUTIONS SDN BHD

DATE:

9/1/19