

KEMENTERIAN KESIHATAN MALAYSIA

MEET Planned Preventive Maintenance Checklist WORK ORDER NO . Micromotor CHECKLISTNO:CL-096-000 REV.000 BE CODE : DE-019 PART 6 QUANTITATIVE TASKS Tick ($\sqrt{\ }$) where appropriate Description Units / Set Measured UOM Limit/Tolerance Values Speed Values PASS FAIL rpm 35,000 ± 5% ()()() PART 7 ELECTRICAL SAFETY TEST ELECTRICAL SAFETY TEST, (attach report) Tick (\forall) where appropriate Standard use : | IEC 60601 | IEC 61010 | IEC 62353 Result: PARISE NOTES FAIL pper done. CORRECTIVE MAINTENANCE REQUIRED FUNCTIONING WORK ORDER NO .__ NOT FUNCTIONING NEXT PPM DATE - 11/2019 PPM has been performed in accordance to the checklist and the equipment is functioning to the intended purpose. COMPLETED BY: FOSTER ANAK ENGKASAN QUANTUM BIOMEDICAL ENGINEERING KUCHING SARAWAK. DATE: 27/4/18